

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 505

Burettes, 50 ml, lat. valve., Schellb., AS, conf. cert.




Batch: 43-20

Serial:

Measuring Nr. 42239

Batch size: 20

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9666 ml	29,9922 ml	10,0192 ml
Deviation:	0,0051 ml	0,0151 ml	0,0130 ml
Precision (CV):	0,0102 %	0,0503 %	0,1303 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	77
Date:	09.11.2020
Balance:	Sartorius - MSE224S-0CE -DU

DE-M