

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 505

Burettes, 50 ml, lat. valve., Schellb., AS, conf. cert.




Batch: 06-21

Serial:

Measuring Nr. 53505

Batch size: 30

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9902 ml	30,0110 ml	9,9991 ml
Deviation:	0,0073 ml	0,0148 ml	0,0166 ml
Precision (CV):	0,0146 %	0,0494 %	0,1657 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	25.03.2022
Balance:	Sartorius - LE244S-OCE

