

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 505

Burettes, 50 ml, lat. valve., Schellb., AS, conf. cert.




Batch: 20-23

Serial:

Measuring Nr. 64373

Batch size: 20

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9860 ml	29,9833 ml	10,0049 ml
Deviation:	0,0168 ml	0,0092 ml	0,0158 ml
Precision (CV):	0,0336 %	0,0307 %	0,1580 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	07.09.2023
Balance:	Sartorius - MSE224S-0CE -DU

