

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 902

**Burette, 100 ml, str. stopcock, Schellb., AS, conf. cert.**




**Batch:** 18-21

**Serial:**

**Measuring Nr.:** 47107

**Batch size:** 10

Volume:	100,000 ml	60,000 ml	20,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	99,9188 ml	59,9761 ml	19,9983 ml
<b>Deviation:</b>	0,0019 ml	0,0027 ml	0,0042 ml
<b>Precision (CV):</b>	0,0019 %	0,0045 %	0,0212 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	07.06.2021
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

