

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 905

**Burette, 100 ml, lat. stopc., Schellbach, AS, conf. cert.**




**Batch:** 18-21

**Serial:**

**Measuring Nr.** 53546

**Batch size:** 10

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	100,0072 ml	60,0066 ml	19,9946 ml
Deviation:	0,0198 ml	0,0390 ml	0,0296 ml
Precision (CV):	0,0198 %	0,0650 %	0,1481 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	30.03.2022
Balance:	Sartorius - LE244S-OCE

