

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 905

Burette, 100 ml, lat. stopc., Schellbach, AS, conf. cert.

Batch: 10-24




Serial:

Measuring Nr. 74549

Batch size:

6

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	99,9788 ml	59,9921 ml	20,0062 ml
Deviation:	0,0161 ml	0,0181 ml	0,0164 ml
Precision (CV):	0,0161 %	0,0302 %	0,0822 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

10.09.2025

Balance:

Sartorius -
MSE224S-0CE
-DU

