

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 102

Burettes, 10 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 19-22

Serial:

Measuring Nr. 58865

Batch size: 10

| Volume: | 10,000 ml | 6,000 ml | 2,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value: | 9,9839 ml | 5,9852 ml | 1,9993 ml |
| Deviation: | 0,0012 ml | 0,0017 ml | 0,0026 ml |
| Precision (CV): | 0,0124 % | 0,0288 % | 0,1284 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|-----------------------------------|
| Operator: | 4 |
| Date: | 15.09.2022 |
| Balance: | Sartorius - MSE224S-0CE -DU |

