

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 103

Burettes, 10 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 05-20

Serial:

Measuring Nr. 36485

Batch size: 19

| Volume: | 10,000 ml | 6,000 ml | 2,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value: | 9,9950 ml | 5,9835 ml | 1,9980 ml |
| Deviation: | 0,0033 ml | 0,0009 ml | 0,0049 ml |
| Precision (CV): | 0,0332 % | 0,0154 % | 0,2436 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|-----------------------------------|
| Operator: | 7 |
| Date: | 17.03.2020 |
| Balance: | Sartorius - MSE224S-0CE -DU |

