

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 103

Burettes, 10 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 11-22

Serial:

Measuring Nr. 56345

Batch size: 20

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9943 ml	5,9917 ml	1,9984 ml
Deviation:	0,0063 ml	0,0044 ml	0,0009 ml
Precision (CV):	0,0631 %	0,0738 %	0,0453 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	31.05.2022
Balance:	Sartorius - LE244S-OCE

