

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 103

Burettes, 10 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 40-25

Serial:

Measuring Nr. 76322

Batch size: 40

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9970 ml	5,9916 ml	1,9956 ml
Deviation:	0,0062 ml	0,0040 ml	0,0012 ml
Precision (CV):	0,0618 %	0,0675 %	0,0617 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

27.01.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M