

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 014 105

**Burette, 10 ml, lat. needle valve, amber, AS, conf. cert.**




**Batch:** 19-22

**Serial:**

**Measuring Nr.** 68943

**Batch size:** 2

Volume:	10,000 ml	6,000 ml	2,000 ml
<b>Tolerance:</b>	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
<b>Average value:</b>	9,9950 ml	5,9986 ml	1,9976 ml
<b>Deviation:</b>	0,0009 ml	0,0019 ml	0,0021 ml
<b>Precision (CV):</b>	0,0085 %	0,0319 %	0,1065 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	09.09.2024
<b>Balance:</b>	Sartorius - LE244S-OCE

