

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 252

Burettes, 25 ml, str. stopc., amber gl., AS, conf. cert.




Batch: 48-24

Serial:

Measuring Nr. 76180

Batch size: 10

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0061 ml	15,0065 ml	5,0008 ml
Deviation:	0,0113 ml	0,0104 ml	0,0068 ml
Precision (CV):	0,0452 %	0,0691 %	0,1351 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	15.01.2026
Balance:	Sartorius - MSE224S-0CE -DU

