

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 253

Burettes, 25 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 29-17

Serial:

Measuring Nr. 14091

Batch size: 100

| Volume: | 25,0000 ml | 15,0000 ml | 5,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0300 ml | ± 0,0300 ml | ± 0,0300 ml |
| Average value: | 24,9904 ml | 14,9870 ml | 4,9954 ml |
| Deviation: | 0,0079 ml | 0,0056 ml | 0,0071 ml |
| Precision (CV): | 0,0316 % | 0,0372 % | 0,1425 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 4 |
| Date: | 13.09.2017 |
| Balance: | Sartorius - LE244S-OCE |

