

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 502 L

Burette, 50 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 03-17

Serial:

Measuring Nr. 29874

Batch size: 5

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0146 ml	29,9683 ml	9,9946 ml
Deviation:	0,0042 ml	0,0034 ml	0,0135 ml
Precision (CV):	0,0083 %	0,0115 %	0,1354 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.05.2019
Balance:	Sartorius - LE244S-OCE

