

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 503

Burettes, 50 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 13-25

Serial:

Measuring Nr. 75148

Batch size: 30

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0175 ml	29,9867 ml	10,0095 ml
Deviation:	0,0074 ml	0,0071 ml	0,0104 ml
Precision (CV):	0,0147 %	0,0236 %	0,1043 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 03.11.2025

Balance: Sartorius -
LE244S-OCE

DE-M