

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 014 504

**Burettes, 50 ml, lat. stopc., amber glass, AS, conf. cert.**




**Batch:** 38-17

**Serial:**

**Measuring Nr.** 15517

**Batch size:** 4

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	49,9924 ml	29,9776 ml	9,9922 ml
<b>Deviation:</b>	0,0044 ml	0,0005 ml	0,0080 ml
<b>Precision (CV):</b>	0,0089 %	0,0017 %	0,0798 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	27.10.2017
<b>Balance:</b>	Sartorius - LE244S-OCE

