

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 505

Burettes, 50 ml, lat. stopc., amber glass, AS, conf. cert.




Batch: 38-17

Serial:

Measuring Nr. 37739

Batch size: 23

| Volume: | 50,000 ml | 30,000 ml | 10,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 50,0368 ml | 29,9920 ml | 10,0000 ml |
| Deviation: | 0,0044 ml | 0,0098 ml | 0,0093 ml |
| Precision (CV): | 0,0087 % | 0,0328 % | 0,0930 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 06.05.2020 |
| Balance: | Sartorius - LE244S-OCE |

