

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 022 057

Burette, 5 ml, straight needle valve, AS, conf. cert.




Batch: 40-18

Serial:

Measuring Nr. 27532

Batch size: 2

| Volume: | 5,0000 ml | 3,0000 ml | 1,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value: | 5,0035 ml | 3,0038 ml | 1,0024 ml |
| Deviation: | 0,0028 ml | 0,0017 ml | 0,0009 ml |
| Precision (CV): | 0,0566 % | 0,0557 % | 0,0936 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 07.02.2019 |
| Balance: | Sartorius - LE244S-OCE |

