

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 022 057

**Burette, 5 ml, straight needle valve, AS, conf. cert.**

**Batch:** 40-18




**Serial:**

**Measuring Nr.** 31630

**Batch size:**

1

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Measured value:	5,0064 ml	3,0043 ml	1,0020 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

16.07.2019

**Balance:**

Sartorius -  
LE244S-OCE

DE-M