

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 022 057

Burette, 5 ml, straight needle valve, AS, conf. cert.




Batch: 12-20

Serial:

Measuring Nr. 44631

Batch size: 2

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0081 ml	3,0027 ml	1,0027 ml
Deviation:	0,0007 ml	0,0011 ml	0,0001 ml
Precision (CV):	0,0142 %	0,0379 %	0,0142 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.02.2021
Balance:	Sartorius - LE244S-OCE

