

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 022 257

Burette, 25 ml, straight needle valve, AS, conf. cert.




Batch: 31-21

Serial:

Measuring Nr.: 49346

Batch size: 49

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0027 ml	15,0196 ml	5,0146 ml
Deviation:	0,0102 ml	0,0030 ml	0,0049 ml
Precision (CV):	0,0407 %	0,0198 %	0,0980 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	28.09.2021
Balance:	Sartorius - LE244S-OCE

