

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 022 507

**Burettes, 50 ml, straight needle valve, AS, conf. cert.**




**Batch:** 13-19

**Serial:**

**Measuring Nr.** 39993

**Batch size:** 2

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0209 ml	30,0181 ml	10,0074 ml
Deviation:	0,0009 ml	0,0005 ml	0,0127 ml
Precision (CV):	0,0017 %	0,0017 %	0,1270 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	09.09.2020
Balance:	Sartorius - LE244S-OCE

