

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 024 107

Burette, 10 ml, str. needle valve, amber, AS, conf. cert.




Batch: 21-16

Serial:

Measuring Nr.: 34655

Batch size: 2

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0081 ml	6,0012 ml	2,0011 ml
Deviation:	0,0044 ml	0,0075 ml	0,0031 ml
Precision (CV):	0,0440 %	0,1254 %	0,1562 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	26.11.2019
Balance:	Sartorius - MSE224S-0CE -DU

