

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 105

Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.




Batch: 20-23

Serial:

Measuring Nr. 63694

Batch size: 51

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0050 ml	5,9924 ml	2,0030 ml
Deviation:	0,0026 ml	0,0026 ml	0,0038 ml
Precision (CV):	0,0261 %	0,0426 %	0,1924 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.07.2023
Balance:	Sartorius - LE244S-OCE

