

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 107

Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.




Batch: 33-17

Serial:

Measuring Nr. 61420

Batch size: 2

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9989 ml	5,9886 ml	1,9962 ml
Deviation:	0,0040 ml	0,0016 ml	0,0074 ml
Precision (CV):	0,0397 %	0,0272 %	0,3690 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	30.01.2023
Balance:	Sartorius - LE244S-OCE

