

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 108

Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.




Batch: 07-20

Serial:

Measuring Nr. 39390

Batch size: 10

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0158 ml	6,0077 ml	1,9987 ml
Deviation:	0,0020 ml	0,0037 ml	0,0030 ml
Precision (CV):	0,0197 %	0,0623 %	0,1519 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	07.08.2020
Balance:	Sartorius - LE244S-OCE

DE-M