

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 108

**Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.**

**Batch:** 02-24




**Serial:**

**Measuring Nr.** 66682

**Batch size:**

1

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Measured value:	9,9996 ml	5,9877 ml	1,9920 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

05.04.2024

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

DE-M