

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 256

Pellet-burettes, 25 ml, lateral stopcock, AS, conf. cert.

Batch: 38-15




Serial:

Measuring Nr.: 19423

Batch size:

6

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9936 ml	14,9986 ml	5,0153 ml
Deviation:	0,0116 ml	0,0117 ml	0,0069 ml
Precision (CV):	0,0463 %	0,0783 %	0,1375 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

04.04.2018

Balance:

Sartorius -
LE244S-OCE

DE-M