

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 504

Pellet-burettes, 50 ml, lateral stopc., AS, conf. cert.

Batch: 18-26




Serial:

Measuring Nr. 78515

Batch size:

10

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	$\pm 0,0500$ ml	$\pm 0,0500$ ml	$\pm 0,0500$ ml
Average value:	50,0176 ml	30,0145 ml	9,9960 ml
Deviation:	0,0023 ml	0,0059 ml	0,0056 ml
Precision (CV):	0,0046 %	0,0196 %	0,0561 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

16.06.2026

Balance:

Sartorius -
LE244S-OCE

DE-M