

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 506

Pellet-burette, 50 ml, lat. stopc., glass, AS, conf. cert.




Batch: 29-20

Serial:

Measuring Nr. 50307

Batch size: 10

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0332 ml	30,0166 ml	10,0064 ml
Deviation:	0,0078 ml	0,0094 ml	0,0104 ml
Precision (CV):	0,0156 %	0,0313 %	0,1042 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	09.11.2021
Balance:	Sartorius - MSE224S-0CE -DU

