

Certificate of performance

Article: 3 122 508

Pellet-burettes, 50 ml, lat. needle valve, AS, conf. cert.




Batch: 29-20

Serial:

Measuring Nr. 56554

Batch size: 20

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0249 ml	30,0105 ml	10,0134 ml
Deviation:	0,0125 ml	0,0054 ml	0,0120 ml
Precision (CV):	0,0250 %	0,0182 %	0,1205 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	13.06.2022
Balance:	Sartorius - MSE224S-0CE -DU