

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 905

Pellet-burettes, 100 ml, lat. needle valve, AS, conf. cert.




Batch: 10-15

Serial:

Measuring Nr. 19187

Batch size: 10

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	100,0059 ml	60,0181 ml	20,0352 ml
Deviation:	0,0413 ml	0,0164 ml	0,0140 ml
Precision (CV):	0,0413 %	0,0273 %	0,0702 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	27.03.2018
Balance:	Sartorius - LE244S-OCE

DE-M