

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 905

**Pellet-burettes, 100 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 10-15

**Serial:**

**Measuring Nr.** 22130

**Batch size:** 4

<b>Volume:</b>	<b>100,000 ml</b>	<b>60,000 ml</b>	<b>20,000 ml</b>
<b>Tolerance:</b>	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
<b>Average value:</b>	100,0033 ml	60,0352 ml	20,0531 ml
<b>Deviation:</b>	0,0356 ml	0,0238 ml	0,0080 ml
<b>Precision (CV):</b>	0,0356 %	0,0397 %	0,0400 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	10.07.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

