

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 105

Pellet-burettes, 10 ml, amber glass, AS, conf. cert.




Batch: 23-21

Serial:

Measuring Nr. 77951

Batch size: 20

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	$\pm 0,0200$ ml	$\pm 0,0200$ ml	$\pm 0,0200$ ml
Average value:	9,9961 ml	6,0045 ml	2,0050 ml
Deviation:	0,0033 ml	0,0043 ml	0,0032 ml
Precision (CV):	0,0332 %	0,0720 %	0,1575 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

05.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M