

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 255

Pellet-burettes, 25 ml, amber glass, AS, conf.cert.

Batch: 22-13




Serial:

Measuring Nr. 21425

Batch size:

7

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9911 ml	14,9936 ml	5,0038 ml
Deviation:	0,0030 ml	0,0060 ml	0,0069 ml
Precision (CV):	0,0119 %	0,0398 %	0,1389 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	04.06.2018
Balance:	Sartorius - LE244S-OCE

DE-M