

Certificate of performance



Article: 3 124 258

Pellet-burettes, 25 ml, amber glass, AS, conf. cert.




Batch: 06-14

Serial:

Measuring Nr.: 21421

Batch size: 8

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9809 ml	14,9942 ml	4,9958 ml
Deviation:	0,0008 ml	0,0066 ml	0,0053 ml
Precision (CV):	0,0031 %	0,0438 %	0,1069 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	04.06.2018
Balance:	Sartorius - LE244S-OCE

