

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 258

**Pellet-burettes, 25 ml, amber glass, AS, conf. cert.**




**Batch:** 48-22

**Serial:**

**Measuring Nr.:** 66155

**Batch size:** 20

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9908 ml	14,9888 ml	5,0076 ml
Deviation:	0,0067 ml	0,0079 ml	0,0070 ml
Precision (CV):	0,0267 %	0,0526 %	0,1399 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

27.02.2024

**Balance:**

Sartorius -  
LE244S-OCE

DE-M