

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 504

**Pellet-burette, 50 ml, lat. stopc., amber, AS, conf. cert.**




**Batch:** 32-18

**Serial:**

**Measuring Nr.:** 25639

**Batch size:** 5

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9905 ml	29,9619 ml	9,9897 ml
Deviation:	0,0078 ml	0,0029 ml	0,0119 ml
Precision (CV):	0,0157 %	0,0097 %	0,1192 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	15.11.2018
Balance:	Sartorius - LE244S-OCE

