

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 504

**Pellet-burette, 50 ml, lat. stopc., amber, AS, conf. cert.**

**Batch:** 43-21




**Serial:**

**Measuring Nr.:** 52772

**Batch size:**

9

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0209 ml	30,0117 ml	10,0188 ml
Deviation:	0,0076 ml	0,0082 ml	0,0123 ml
Precision (CV):	0,0152 %	0,0273 %	0,1229 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

03.03.2022

**Balance:**

Sartorius -  
LE244S-OCE

DE-M