

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 505

Pellet-burettes, 50 ml, amber glass, AS, conf. cert.




Batch: 29-19

Serial:

Measuring Nr. 34079

Batch size: 24

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9753 ml	29,9851 ml	10,0149 ml
Deviation:	0,0085 ml	0,0145 ml	0,0122 ml
Precision (CV):	0,0169 %	0,0485 %	0,1219 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	04.11.2019
Balance:	Sartorius - LE244S-OCE

