

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 508

**Pellet-burettes, 50 ml, amber glass, AS, conf. cert.**




**Batch:** 25-20

**Serial:**

**Measuring Nr.:** 42578

**Batch size:** 10

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0007 ml	29,9854 ml	10,0116 ml
Deviation:	0,0173 ml	0,0121 ml	0,0032 ml
Precision (CV):	0,0346 %	0,0403 %	0,0318 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	16.11.2020
Balance:	Sartorius - MSE224S-0CE -DU

