

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 152 505

Univ. autom. burette, 50 ml, lat. valve, AS, conf. cert.




Batch: 29-20

Serial:

Measuring Nr. 66152

Batch size: 2

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9724 ml	29,9915 ml	9,9927 ml
Deviation:	0,0094 ml	0,0027 ml	0,0116 ml
Precision (CV):	0,0187 %	0,0090 %	0,1159 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	26.02.2024
Balance:	Sartorius - LE244S-OCE

