

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 154 504

**Univ. autom. burette, 50 ml, amber glass, AS, conf. cert**

**Batch:** 29-19




**Serial:**

**Measuring Nr.** 34045

**Batch size:**

6

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0223 ml	29,9840 ml	10,0128 ml
Deviation:	0,0105 ml	0,0169 ml	0,0063 ml
Precision (CV):	0,0211 %	0,0563 %	0,0635 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

31.10.2019

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

