

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 154 504

Univ. autom. burette, 50 ml, amber glass, AS, conf. cert




Batch: 43-21

Serial:

Measuring Nr. 54565

Batch size: 23

| Volume: | 50,0000 ml | 30,0000 ml | 10,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 50,0041 ml | 30,0052 ml | 10,0103 ml |
| Deviation: | 0,0192 ml | 0,0149 ml | 0,0132 ml |
| Precision (CV): | 0,0384 % | 0,0495 % | 0,1322 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

26.04.2022

Balance:

Sartorius -
LE244S-OCE

DE-M