

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 372 012

**Mikro burettes, 1 ml, str. stopc., glass, AS, conf. cert.**




**Batch:** 13-18

**Serial:**

**Measuring Nr.** 20735

**Batch size:** 3

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0051 ml	0,6027 ml	0,1999 ml
Deviation:	0,0003 ml	0,0005 ml	0,0010 ml
Precision (CV):	0,0266 %	0,0754 %	0,4874 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	16.05.2018
Balance:	Sartorius - LE244S-OCE

