

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 012

Mikro burette, 1 ml, str. stopc., glass, AS, conf. cert.




Batch: 13-18

Serial:

Measuring Nr. 66724

Batch size: 9

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0015 ml	0,6030 ml	0,2009 ml
Deviation:	0,0008 ml	0,0013 ml	0,0004 ml
Precision (CV):	0,0811 %	0,2213 %	0,1900 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	05.04.2024
Balance:	Sartorius - LE244S-OCE

