

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 372 022

**Micro burette, 2 ml, str. stopc., glass, AS, conf. cert.**

**Batch:** 41-20




**Serial:**

**Measuring Nr.:** 41821

**Batch size:**

5

| Volume:         | 2,0000 ml   | 1,2000 ml   | 0,4000 ml   |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value:  | 2,0015 ml   | 1,1981 ml   | 0,4005 ml   |
| Deviation:      | 0,0028 ml   | 0,0006 ml   | 0,0021 ml   |
| Precision (CV): | 0,1386 %    | 0,0483 %    | 0,5323 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

28.10.2020

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

