

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 022

Micro burette, 2 ml, str. stopc., glass, AS, conf. cert.




Batch: 13-23

Serial:

Measuring Nr. 66736

Batch size: 10

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	$\pm 0,0100$ ml	$\pm 0,0100$ ml	$\pm 0,0100$ ml
Average value:	1,9997 ml	1,2014 ml	0,4007 ml
Deviation:	0,0034 ml	0,0020 ml	0,0006 ml
Precision (CV):	0,1706 %	0,1674 %	0,1382 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

11.04.2024

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M