

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 023

Micro burette, 2 ml, str. stopc., PTFE, AS, conf. cert.

Batch: 41-20




Serial:

Measuring Nr. 43548

Batch size:

6

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9974 ml	1,1973 ml	0,3991 ml
Deviation:	0,0019 ml	0,0029 ml	0,0017 ml
Precision (CV):	0,0949 %	0,2408 %	0,4156 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	13.01.2021
Balance:	Sartorius - MSE224S-0CE -DU

