

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 052

Micro burettes, 5 ml, straight stopc., AS, conf. cert.




Batch: 12-23

Serial:

Measuring Nr. 67066

Batch size: 10

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0037 ml	3,0057 ml	1,0037 ml
Deviation:	0,0018 ml	0,0017 ml	0,0013 ml
Precision (CV):	0,0352 %	0,0581 %	0,1280 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	02.05.2024
Balance:	Sartorius - LE244S-OCE

DE-M